

HeLP Challenge

- Cardiac Segmentation

engineering advices

Nov. 28. 2018

Jieun Kim

Radiographic presentations

Anterior


Axial plane



Right

Left

Posterior



Head

Coronal plane



Right

Left

Foot

Head

Sagittal plane



Anterior

Posterior

Foot

MHA analyzed file format 1/3

- MHA file
 - An ITK(Insight segmentation and Registration Toolkit) MetaImage
 - Soft application of viewer : ITK SNAP, MITK, ImageJ etc.
- Load MHA file in python

- Using simpleITK

```
import SimpleITK as sitk
```

```
imageInput = sitk.ReadImage(file_name_imageT1[n])
```

- Get size, origin, spacing, direction

```
imageInput = sitk.ReadImage(file_name_imageT1[0])
```

```
origin3d = imageInput.GetOrigin()
```

```
spacing3d = imageInput.GetSpacing()
```

```
size3d = imageInput.GetSize()
```

```
direction3d = imageInput.GetDirection()
```

MHA analyzed file format 2/3

- Convert itk image to array

```
imageT = sitk.GetArrayFromImage(image) # covert itk image to Array  
model.train_on_batch(imageT, LabelT1)
```


- Write image

```
model1.load_weights(path_model1)  
result11 = model1.predict(imageT)  
result1 = result11.reshape((img_rows, img_cols, img_dep))  
ResultImg11 = sitk.GetImageFromArray(result1)  
  
final_img1.SetSpacing(original_img.GetSpacing())  
final_img1.SetOrigin(original_img.GetOrigin())  
  
sitk.WriteImage(final_img1, output_filename)
```

MHA analyzed file format 3/3





- Soft application of visualization

dragging




Pre-processing -intensity windowing

- (WindowLevel, WindowWidth) = (0, 600) ~ (0, 1000)



Reference method flow chart



Reference Method 1/2

```
for epoch in range(20) :
    random.shuffle(indices_t)
    random.shuffle(indices_v)
    for n in indices_t:

        input_img = sitk.ReadImage(file_name_imageT1[n])# read image
        input_lbl = sitk.ReadImage(file_name_imageT2[n])# read label

        isocubic_img = resample_ISOCUBIC(input_img,(0.5,0.5,0.5))
        isocubic_label=resampleMask_ISOCUBIC(input_lbl,(0.5,0.5,0.5)) → sitk.Resample()

        isocubic_img_ww = sitk.Cast(sitk.IntensityWindowing(isocubic_img, LowT,HighT,0,255), sitk.sitkUInt8)

        size3d = isocubic_img_ww.GetSize()
        size2d = (size3d[0],size3d[1],0) #axial slices,

        outimg = sitk.Image(isocubic_img.GetSize(), sitk.sitkUInt8)
        outimg.SetDirection(isocubic_img.GetDirection())
        outimg.SetSpacing(isocubic_img.GetSpacing())
        outimg.SetOrigin(isocubic_img.GetOrigin())

        outimage_array = sitk.GetArrayFromImage(outimg)

        num_of_slices = isocubic_img_ww.GetDepth()
```


Reference Method 2/2

```
for i in range(num_of_slices):
    prev_idx = (0,0,i-displacement)
    cur_idx = (0,0,i)
    next_idx = (0,0,i+displacement)

    if( prev_idx[2] <0 or (next_idx[2] >=num_of_slices):
        continue

    pre_image = sitk.Extract(isocubic_img_ww,size2d,prev_idx)
    cur_image = sitk.Extract(isocubic_img_ww,size2d,cur_idx)
    next_image = sitk.Extract(isocubic_img_ww,size2d,next_idx)

    mask2d = sitk.Extract(isocubic_label,size2d, cur_idx)

    #512*512로 추출한 영상들 resize하기.
    prev_image_resample = resample_resize2D(prev_image,(img_rows,img_cols))
    cur_image_resample = resample_resize2D(cur_image,(img_rows,img_cols))
    next_image_resample = resample_resize2D(next_image,(img_rows,img_cols))

    mask2d_resample = resample_resizeMask2D(mask2d,(img_rows,img_cols)) # training할 최종 label
    mask_final = label_Mask_3ch(mask2d_resample) # 3label을 각각 label별로 채널에 나누어서 저장.
    rgbImage = sitk.Compose(prev_image_resample, cur_image_resample,next_image_resample)


    image_train_arr = sitk.GetArrayFromImage(rgbImage)
    label_train_arr = sitk.GetArrayFromImage(mask_final)

    image_training = image_train_arr.reshape((1,img_rows,img_cols,3) )
    label_training = label_train_arr.reshape((1,img_rows,img_cols,n_class))


    model.train_on_batch( image_training, label_training)
```

Result of ref. method


- Myocardium segmentation in cardiac CT images




Manual Segmentation



Multi-atlas Segmentation



FCN Segmentation



$$DSC(A, B) = \frac{2 * |A \cap B|}{|A| + |B|}$$

Dice Coefficient result	Min	Max	Mean	Std
Multi-atlas	74.36674	92.10817	87.73908	3.058546
FCN	84.80109	95.0653	90.84953	2.264284

Running Time

► Multi-atlas : 40min/1case

► FCN : 2~3min/1case


TIPS : Docker image build (nvidia-docker)

- If your cuda is less than 9.0 version, do not use `FROM tensorflow/tensorflow:latest-gpu`
 - It requires **CUDA 9.0 or more.** (Nov.27. 2018)
 - 1.2.0-gpu-py3
- Make good use of official docker-hub
 - <https://www.tensorflow.org/install/docker>
 - <https://hub.docker.com/r/nvidia/cuda/>
- Always commit docker-image when you modify your container.

Command : Docker image build (nvidia-docker)

- `nvidia-docker exec -it [container-id] bash`
- `docker kill [container id] # remove container`
- `docker kill $(docker ps -q) #remove all container`
- `docker rmi -f [image id] # remove an image`
- `docker system prune -a # remove all image`
- `docker cp [host file or folder name] [container id]:[container file or folder name]`
copy files
- `docker commit [container name] [img name] # commit!`
- `docker save [img name]:[img tag] | qzip > [filename] # docker save image (.tar.gz)`
- `qunzip -c [filename] | docker load # docker qzip file load`

Dataset structure



- [train]

- /data/train/CHD/image
- /data/train/CHD/label
- /data/train/HCMP/image
- /data/train/HCMP/label

- [test]

- - /data/test/CHD
- - /data/test/HCMP

- [output]

- /data/output/CHD
 - '%03d_dia_output.mha'
 - '%03d_sys_output.mha'
- /data/output/HCMP
 - 'A_%03d_output.mha'
 - 'N_%03d_output.mha'



Outline

-cardiac segmentation




- [output]
 - 3-dimensional-1channel MHA file, (x,y,z,1)
 - Having the same structure as the test image
 - Size, origin, spacing, direction
 - Each output consisting of 3 labels
- To evaluate the performance, Dice similarity coefficient (DSC) of the cardiac segmentation will be computed.
- $DSC(R, G) = 2 |R \cap G| / (|R| + |G|)$

Cardiac dataset (phase1)

- Three-dimensional CT scans for systolic, diastolic, abnormal and normal, respectively.
 - Train : 70 cases, test : 10cases (diastolic and systolic, respectively) → CHD
 - Train : 50 cases, test : 7 cases (abnormal)  HCMP
 - Train : 12 cases, test : 3 cases (normal)  HCMP
- The image resolution is 512×512 voxels in-plane and the range of all slice number is diverse.
- Various range of voxel spacing
- The train dataset includes the segmentation result by experts to label and delineate

Cardiac CT data 1/2 -CHD

Diastolic




Systolic






- LVM = left ventricle myocardium(1, red)
- LV = left ventricle chamber (2, green)
- RV = right ventricle chamber (3, blue)



Diastolic






Diastolic






Systolic



Systolic





Short axis view - diastolic





Cardiac CT data 2/2 -HCMP

Normal






Abnormal






- LVM = left ventricle myocardium (1, red)
- APM = anterior papillary muscle (2, green)
- PPM = posteromedial papillary muscle (3, blue)


Normal



Normal



Abnormal



Abnormal

